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1. Report of japan neurosurgery registry (2015 - 2017). Neurologia medico-chirurgica. 2019;59(Special-Issue):13-81.10.2176/nmc.si.2019-0001
2. Brott TG, Halperin JL, Abbara S, et al.: 2011 asa/accf/aha/aann/aans/acr/asnr/cns/saip/scail/sir/snis/svm/svs guideline on the management of patients with extracranial carotid and vertebral artery disease: Executive summary: A report of the american college of cardiology foundation/ american heart association task force on practice guidelines, and the american stroke association, american association of neuroscience nurses, american association of neurological surgeons, american college of radiology, american society of neuroradiology, congress of neurological surgeons, society of atherosclerosis imaging and prevention, society for cardiovascular angiography and interventions, society of interventional radiology, society of neurointerventional surgery, society for vascular medicine, and society for vascular surgery. J Am Coll Cardiol. 2011;57(8):1002-1044.10.1016/j.jacc.2010.11.005
3. AbuRahma AF, Robinson PA, Saiedy S, et al.: Prospective randomized trial of carotid endarterectomy with primary closure and patch angioplasty with saphenous vein, jugular vein, and polytetrafluoroethylene: Long-term follow-up. J Vasc Surg. 1998;27(2):222-232; discussion 233-224.10.1016/s0741-5214(98)70353-2
4. Hertzler NR, O'Hara PJ, Mascha EJ, et al.: Early outcome assessment for 2228 consecutive carotid endarterectomy procedures: The cleveland clinic experience from 1989 to 1995. J Vasc Surg. 1997;26(1):1-10.10.1016/s0741-5214(97)70139-3
5. Orrapin S, Benyakorn T, Howard DP, et al.: Patches of different types for carotid patch angioplasty. Cochrane Database Syst Rev. 2021;2(2):CD000071.10.1002/14651858.CD000071.pub4
6. 遠藤 俊, 飯原 弘, 永田 泉, et al.: 頸動脈高度狭窄病変の本邦治療の現状 : j c a s から. 脳卒中. 2005;27(4):492-497.10.3995/jstroke.27.492
7. Huizing E, Vos CG, van den Akker PJ, et al.: A systematic review of patch angioplasty versus primary closure for carotid endarterectomy. J Vasc Surg. 2019;69(6):1962-1974.e1964.10.1016/j.jvs.2018.10.096

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